

SKIN FOR DIRECT IMMUNOFLUORESCENCE

Indications:

For detection of disease in skin biopsies.

Specimen:

Skin biopsies

Collection:

- Send biopsy in Michel's media. Michel's media can be requested prior to the procedure. If Michel's media is not available at the time of biopsy, send biopsy on saline-moistened gauze, and the specimen will be transferred to Michel's media at the time of receipt.
 - Note: Call courier dispatch at (702) 249-4094. ARDX LMC can deliver Michel's media when called 24 hours prior to the biopsy.
- Label specimen container with patient name, specimen type and location.
- Submit all patient information following the procedure for "Completing a Surgical Requisition" (page 5).