

## **PROSTATE SEXTANT BIOPSY KITS**

Urology prostate sextant biopsy kits will be provided for clients. This kit will include 8 containers of formalin for needle biopsies from eight different sites: Rt. Base, Rt. Mid, Rt. Apex, Rt. Lat Mid, Lt. Base, Lt. Mid, Lt. Apex, and Lt. Lat. Mid.

Each specimen container will have a color-coded label on the lid stating the site:

RIGHT BASE	RED
RIGHT MID	PINK
RIGHT APEX	ORANGE
RIGHT LATERAL MID	PURPLE
LEFT BASE	YELLOW
LEFT MID	BLUE
LEFT APEX	GREEN
LEFT LATERAL MID	DARK BLUE

The client will need to label the container on the side with the patient name and place in the appropriate location holder in the kit. Unused vials should not be returned. The requisition should be folded and placed in the plastic pouch on the lid of the box.